



Weill Cornell Medicine

ApHID: pH-Sensitive Fluorescent Dye Resistant to Photobleaching and Oxidation

Lead Inventors:

David Warren, Ph.D.

Adjunct Associate Professor of Research in
Biochemistry, Biochemistry, Weill Cornell Medicine

Frederick Maxfield, Ph.D.

Vladimir Horowitz and Wanda Toscanini Horowitz
Distinguished Professor in Neuroscience,
Biochemistry, Weill Cornell Medicine
Professor of Biochemistry, Biochemistry, Weill
Cornell Medicine

Santiago Sole Domenech, Ph.D.

Assistant Professor of Research in Biochemistry,
Biochemistry, Weill Cornell Medicine

Business Development Contact:

Jamie Brisbois

Manager, Business Development and Licensing

(646) 962-7049

jamie.brisbois@cornell.edu

ApHID: pH-Sensitive Fluorescent Dye Resistant to Photobleaching and Oxidation

Background & Unmet Need

- Human intercellular activities are often facilitated by movement of membrane-bound vesicles such as late endosomes and lysosomes (LE/Ly)
- LE/Ly dysfunction has been linked to several diseases and disorders including Alzheimer's disease and Tay-Sachs disease
- The ability to precisely measure LE/Ly activities under different cellular environments is important for the understanding of these diseases, and their potential diagnosis
- However, existing commercial fluorescence lack the photo intensity and chemical stability in highly acidic and reactive cellular environments within LE/Ly vesicles
- **Unmet Need:** A fluorescence probe that is pH-sensitive and able to withstand oxidation and photobleaching while maintaining structural integrity *in vivo*

Technology Overview

- **The Technology:** Acidic pH indicator Dye (ApHID) with high resistance to oxidation and photobleaching
- ApHID is composed of a BODIPY core and determines pH of the environment using an aniline moiety that has two methyl groups attached
- Optimized for use between pH 4.0 – 6.0, ApHID's fluorescence emission increases sharply in amplitude with increasing acidity
- ApHID has pKa of 5.4 and excitation max at 506 nm
- **PoC Data:** ApHID fluorescence is 12-fold greater at pH 4.0 relative to pH 6.0
- ApHID fluorescence output showed no decrease after photobleaching in live cells and decreased by only 12% in fixed cells, compared to an 83% and 82% decrease with fluorescein and Oregon Green
- ApHID exhibits the greatest fluorescent dynamic range at the physiological pH range of LE/Lys compared to currently available commercial dyes

Inventors:

David Warren
Frederick Maxfield
Santiago Sole Domenech

Patents:

[PCT Application Filed](#)

Publications:

[Solé-Domènech et al. Cell Reports Methods. 2025.](#)

Biz Dev Contact:

Jamie Brisbois
(646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9361



Weill Cornell Medicine

ApHID: pH-Sensitive Fluorescent Dye Resistant to Photobleaching and Oxidation

Technology Applications

- Fluorescent dye for LE/Ly research and experiments, including pH and live cell imaging
- Labeling of A β aggregates to image digestive exophagy
- Tracking efficacy of drugs for neurodegenerative diseases
- Tool for cancer research and drug development

Technology Advantages

- Greater brightness and sensitivity to acidity than existing dyes
- Highly resistant to photobleaching and highly concentrated reactive oxygen species
- Stable in living cells over 15 hours+ while emitting strong fluorescent signal
- Enables detection of pH differences within cells and real-time changes in pH

Supporting Data / Figures




Figure 1: (A) Chemical structure of ApHID. (B) Excitation and (C) emission spectra plotted against buffer pH, measured for a 0.04 mg/mL dilution of 10 KDa amino-dextran labeled with ApHID at a 1.6:1 molar ratio, in pH-adjusted buffers.

Inventors:

David Warren
Frederick Maxfield
Santiago Sole Domenech

Patents:

[PCT Application Filed](#)

Publications:

[Solé-Domènech et al. Cell Reports Methods. 2025.](#)

Biz Dev Contact:

Jamie Brisbois
(646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9361



Weill Cornell Medicine

ApHID: pH-Sensitive Fluorescent Dye Resistant to Photobleaching and Oxidation

Supporting Data / Figures




Figure 2: (A, B, C) ApHID, fluorescein, Oregon Green, and LysoSensor y/b and/or BioTracker, pHrodo absorbance and emission profiles vs buffer pH. (D) ApHID, Oregon Green, and fluorescein fluorescence titrations against buffer pH in the presence of hydroxyl radical. (E-F) ApHID, Oregon Green, and fluorescein fluorescence in high salt environment.

Inventors:

David Warren
Frederick Maxfield
Santiago Sole Domenech

Patents:

[PCT Application Filed](#)

Publications:

[Solé-Domènech et al. Cell Reports Methods. 2025.](#)

Biz Dev Contact:

Jamie Brisbois
(646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9361



Weill Cornell Medicine

ApHID: pH-Sensitive Fluorescent Dye Resistant to Photobleaching and Oxidation

Supporting Data / Figures




Figure 3: (A-F) Comparative photostability of ApHID, fluorescein, and Oregon Green LE/Lys of fixed J774 macrophages, imaged by confocal microscopy (G-J) Comparative photostability of ApHID and LysoSensor yellow/blue in live J774 macrophages (K) Normalized fluorescence intensity plotted against irradiation cycle for ApHID, fluorescein, and Oregon Green dextran in fixed cells and ApHID and LysoSensor in live cells.

Inventors:

David Warren
Frederick Maxfield
Santiago Sole Domenech

Patents:

[PCT Application Filed](#)

Publications:

[Solé-Domènech et al. *Cell Reports Methods*. 2025.](#)

Biz Dev Contact:

Jamie Brisbois
(646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9361



Weill Cornell Medicine

ApHID: pH-Sensitive Fluorescent Dye Resistant to Photobleaching and Oxidation

Supporting Data / Figures





Figure 4: Left: (J-M) ApHID enabled detection of small differences in pH between individual compartments within single cells (N-O) ApHID is sensitive to changes in pH in real time. A brief treatment with 20 mM methylamine induced rapid alkalinization **Right:** ApHID reports LE/Ly pH in a stable manner over time.

Inventors:

David Warren
Frederick Maxfield
Santiago Sole Domenech

Patents:

[PCT Application Filed](#)

Publications:

[Solé-Domènech et al. *Cell Reports Methods*. 2025.](#)

Biz Dev Contact:

Jamie Brisbois
(646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9361

ApHID: pH-Sensitive Fluorescent Dye Resistant to Photobleaching and Oxidation

Supporting Data / Figures






Figure 5: Intravital LE/Ly pH imaging of macrophages using ApHID allows for detection of dynamic changes in organelle acidification in real time. Pro-inflammatory ischemic lesions were induced in the parenchyma of mice via photothrombosis and LE/Ly acidification in local macrophages was measured. Lesioned areas displayed phagocytic cells with amoeboid morphologies undergoing LE/Ly acidification (red arrowheads). Re-acidification reached its peak 7 days after ischemic insult and seemed to plateau at 11 days post lesion.

Inventors:

David Warren
Frederick Maxfield
Santiago Sole Domenech

Patents:

[PCT Application Filed](#)

Publications:

[Solé-Domènech et al. Cell Reports Methods. 2025.](#)

Biz Dev Contact:

Jamie Brisbois
(646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9361

ApHID: pH-Sensitive Fluorescent Dye Resistant to Photobleaching and Oxidation

Supporting Data / Figures




Figure 6: Left: In vivo 3-photon visualization of LE/Lys using Aphid **Top Right:** ApHID/Alexa 405 fluorescence ratio measured in fixed macrophages, plotted against buffer pH **Bottom Right:** ApHID can be used to monitor acidification, and alkalinization, of LE/Ly pH in vivo and in real time accurately. Mice were treated with chloroquine (CQ), a weak base. Macrophages with acidic LE/Ly pH underwent rapid alkalinization 2 min and 10 min, reported by a sharp drop in ApHID fluorescence.

Inventors:

David Warren
Frederick Maxfield
Santiago Sole Domenech

Patents:

[PCT Application Filed](#)

Publications:

[Solé-Domènech et al. *Cell Reports Methods*. 2025.](#)

Biz Dev Contact:

Jamie Brisbois
(646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9361



Weill Cornell Medicine



Weill Cornell Medicine