



Weill Cornell Medicine

MRD-EDGE: Ultra-Sensitive Detection of Circulating Tumor DNA for Cancer Screening and Diagnosis

Lead Inventors:

Dan Landau, M.D., Ph.D.

Associate Professor of Medicine, Division of Hematology and Medical Oncology
Associate Professor of Physiology and Biophysics,
Weill Cornell Medicine
Core Member, New York Genome Center

Adam Widman, M.D.

Assistant Attending Physician
Memorial Sloan Kettering Cancer Center

Business Development Contact:

Jamie Brisbois
Manager, Business Development and Licensing

(646) 962-7049
jamie.brisbois@cornell.edu

MRD-EDGE: Ultra-Sensitive Detection of Circulating Tumor DNA for Cancer Screening and Diagnosis

Background & Unmet Need

- Liquid biopsy is an emerging noninvasive method for cancer diagnosis and monitoring which involves sequencing blood plasma cell-free DNA (cfDNA) to identify circulating tumor DNA (ctDNA)
- ctDNA detection is of particular interest for evaluating minimal residual disease (MRD), which indicates the lingering presence of cancerous cells after an initial cancer treatment
- Current ctDNA detection methods have inadequate sensitivity in low volume cancer due to the sparsity of ctDNA in blood and usually require a matched tumor sample, which may not be feasible in many clinical settings
- **Unmet Need:** A sensitive noninvasive liquid biopsy platform to accurately detect residual tumor in blood samples at low tumor burden in the tumor-informed or tumor-naïve context

Technology Overview

- **The Technology:** MRD-EDGE is an ultra-sensitive machine learning-guided ctDNA analysis platform for MRD detection in low tumor fraction cancers
- MRD-EDGE incorporates simultaneous profiling of single nucleotide variants (SNV) and copy number variants (CNV) to enhance ctDNA detection
- The deep learning SNV classifier integrates properties of somatic mutations to distinguish ctDNA from sequencing error, enabling ctDNA detection even in the parts per million range and below
- The CNV module couples read-depth denoising with fragmentomics and an allelic imbalance classifier to detect ctDNA even at low aneuploidy levels
- **PoC Data:** MRD-EDGE enabled tracking tumor burden changes in response to immunotherapy in non-small cell lung cancer (NSCLC), ctDNA shedding in precancerous colorectal adenomas, and de novo mutation calling in melanoma, yielding clinically informative tumor fraction monitoring

Inventors:

Dan Landau
Adam Widman
Minita Shah
Cole Khamnei
Jacob Bass

Patents:

[PCT Application Filed](#)
[PCT Application Filed](#)

Publications:

[Widman et al. bioRxiv](#)
[Preprint]. 2022.

Biz Dev Contact:

Jamie Brisbois
(646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9641
D-10093

MRD-EDGE: Ultra-Sensitive Detection of Circulating Tumor DNA for Cancer Screening and Diagnosis

Technology Applications

- Ultrasensitive MRD detection following surgical resection of cancer
- Noninvasive liquid biopsy for cancer screening
- Real-time serial monitoring of therapy response to inform therapeutic optimization
- Patient monitoring during remission for early detection of relapse
- Applicability in a wide range of solid tumors

Technology Advantages

- Ultra-sensitive SNV and CNV detection due to advanced error suppression and radical amplification of ctDNA signal
- ctDNA detection in tumor-informed or tumor-naïve context (without matched tumor tissue)
- Simple Whole Genome Sequencing (WGS) workflow does not require custom panel creation or molecular barcodes and can work with limited input material

Supporting Data / Figures




Figure 1: ctDNA analysis workflow with MRD-EDGE.

Inventors:

Dan Landau
Adam Widman
Minita Shah
Cole Khamnei
Jacob Bass

Patents:

[PCT Application Filed](#)
[PCT Application Filed](#)

Publications:

[Widman et al. bioRxiv \[Preprint\]. 2022.](#)


Biz Dev Contact:

Jamie Brisbois
(646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9641
D-10093

MRD-EDGE: Ultra-Sensitive Detection of Circulating Tumor DNA for Cancer Screening and Diagnosis



Inventors:

Dan Landau
 Adam Widman
 Minita Shah
 Cole Khamnei
 Jacob Bass

Patents:

[PCT Application Filed](#)
[PCT Application Filed](#)

Publications:

[Widman et al. bioRxiv](#)
 [Preprint]. 2022.


Biz Dev Contact:

Jamie Brisbois
 (646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9641
 D-10093

MRD-EDGE: Ultra-Sensitive Detection of Circulating Tumor DNA for Cancer Screening and Diagnosis



Inventors:

Dan Landau
Adam Widman
Minita Shah
Cole Khamnei
Jacob Bass

Patents:

[PCT Application Filed](#)
[PCT Application Filed](#)

Publications:

[Widman et al. bioRxiv](#)
[Preprint]. 2022.

Biz Dev Contact:

Jamie Brisbois
(646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9641
D-10093

MRD-EDGE: Ultra-Sensitive Detection of Circulating Tumor DNA for Cancer Screening and Diagnosis

Supporting Data / Figures





Figure 4. (A) Tumor-informed post-operative ctDNA detection in early-stage colorectal cancer is associated with shorter progression-free survival. **(B)** Serial tumor-naïve ctDNA monitoring with MRD-EDGE reflects changes seen on imaging.

Inventors:

Dan Landau
Adam Widman
Minita Shah
Cole Khamnei
Jacob Bass

Patents:

[PCT Application Filed](#)
[PCT Application Filed](#)

Publications:

[Widman et al. *bioRxiv*](#)
[Preprint]. 2022.

Biz Dev Contact:

Jamie Brisbois
(646) 962-7049
jamie.brisbois@cornell.edu

Cornell Reference:

D-9641
D-10093



Weill Cornell Medicine